



# Enrichment Sample Homogenizer (ESH)

## Intended Use:

The Enrichment Sample Homogenizer is used for automatic homogenization of food enrichment samples with enrichment media.

## Features of the ESH instrument:

- Multi-function digital display/control panel
- Variable speed (4, 6, 8, and 10 strokes)
- Variable time (1 second to 1 hour or open-running/countdown)
- Close-tight window door
- Adjustable blending power
- Side-by-side paddle stop
- Silent, brushless motor
- Security drip tray
- Removable autoclavable paddles

Category	Specification
Useful volume	50–400 mL
Power source	100V~ to 240V~ 50 to 60 Hz
Size (width x depth x height)	26.5 x 42 x 26 cm 10.433 x 16.45 x 10.24 in.
Weight	16 kg
Blending time	1 second—1 hour, and open running / countdown
Blending speed	4, 6, 8, and 10 strokes

[www.merckmillipore.com/biomonitoring](http://www.merckmillipore.com/biomonitoring)

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

The vibrant M, Merck, DiluCult, and Millipore are trademarks of Merck. All other trademarks belonging to third parties are the properties of their respective owners. Copyright © 2017 Merck KGaA. All Rights Reserved. Lit. No. DS5766ENEU | 2017 - 02830 | 03/2017

## Operator and Equipment Safety

- Place the ESH instrument upright on a clean, flat, stable, horizontal surface, away from any source of excessive heat and close to an easily accessible, properly grounded power supply outlet. Do not run the instrument in any other position.
- Do not open the housing without first unplugging the instrument from the power source.
- Do not expose the ESH instrument to liquid. If this happens, immediately switch off and disconnect the pump from the power outlet, and then decontaminate the ESH instrument.
- Before cleaning, shut down the ESH and disconnect it from the power source.
- The power supply must be protected by a fuse below the main connection.

For further details, see the Enrichment Sample Homogenizer User Guide.

## Related Products:

Instrument	Catalog Number
DiluCult™ dilutor instrument	5 42753 0001
DiluCult™ 2 dilutor instrument	5 42760 0001

Find contact information for your country at: [www.merckmillipore.com/offices](http://www.merckmillipore.com/offices)

For Technical Service, please visit: [www.merckmillipore.com/techservice](http://www.merckmillipore.com/techservice)